


BURGUNDY RED



0° - 113kg
6° - 80kg

UL8-WSLA

FEATURES

- Powerchair pulls apart for storage and transport
- Comfortable mid back contoured seat
- Height / angle adjustable seat
- Lift up & adjustable width / height arm rests
- Fold up adjustable height / depth footrest
- Mid wheel drive
- Finger control operation
- High / Low speed settings
- Battery power indicator
- Horn
- Ideal for indoor and outdoor use on even ground
- Safe to transport on planes

SPECIFICATIONS

Overall: 725mm (L) x 545mm (W) x 1180mm (H)
 Turning Circle: 755mm
 Speed: 6km/h
 Maximum range: Up to 12km
 Safe climb angle: 6°
 Maximum climb angle: 6°
 Motor: 2 x 0.5Hp high torque motors (2 Pole)
 Batteries: 2 x 12 Amp minimum (supplied)
 Charger: 2 Amp
 Tyres: Front: Ø155mm / Rear: Ø200mm
 Maximum User Weight: 113kg - on flat ground

Product Weight: 48.1kg - with batteries
 Seat with Controller Weight: 16.5kg
 Front Chassis with Battery Weight: 18.4kg
 Rear Chassis with Drive Train Weight: 13.2kg