

Emperor

BOSS SCOOTER



The BOSS of all scooters. Engineered from the ground up, this premium bariatric scooter boasts an incredible 250kg load capacity. The OT designed bariatric seat is bench style with easy access and couples with multiple driver comfort settings. Added to its amazing performance ability, the Emperor is arguably Australia's No1 choice!



TopGun
Mobility
WHEELS FOR LIFE

Emperor Features

- Industry leading 5 year warranty
- World's best componentry
- Max load capacity 250kg
- Extra wide bariatric bench seat as standard
- Extra wide tyres for driver comfort
- 22.5cm ground clearance
- All around adjustable h/duty suspension
- Safety hand brake complimenting electric brakes
- Safe cornering steering sensors as standard
- High visibility head and taillights
- Adjustable tiller for driver comfort
- Large rear lockable box
- USB port for phone charging – Peace of mind
- Power save mode

Specifications

Motor Power	1300W
Max. Speed	10km/h
Operation Range	Up to 70km
Max. Climbing Grade	10° (depending on weight of user)
Turning Radius	260cm
Ground Clearance	22.5cm (+/- 2.5cm)
Max. Loading Weight	250kg
Overall Dimensions	165 x 84 x 140cm
Wheel Size - Front	15" pneumatic
Wheel Size - Rear	15" pneumatic
Weight	150kg
Battery	12V 70Ah x 2 (lead-acid)
Controller	PG 200A
Brake System	Electromagnetic
Free Wheel Mode	Yes
Suspension	Yes
Rear Mirrors	Left & Right
Charger Size	24V 8A



NOTE: All specifications subject to change without notice.

- The information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery Condition and tyre pressure.