




Washable
Arm & Headrest
Covers


NEW PRODUCT!




POWERED
Headrest and
Lumbar Support



GERMAN ACTUATORS
Lift You
To Your Feet


COMPLIANT TO AS/NZS 4688



For more information
about this product



Foot Rest
Extension

ROSA

4 MOTOR LIFT CHAIR

LC0801 (Shadow Grey)
LC0802 (Ocean Blue)
LC0803 (Ice Pearl)

FEATURES

- Four motor lift chair
- Lift / Rest / Recline armchair
- Seat raise and tilt forward function
- Backrest recline function
- Legrest elevation function
- Fabric upholstery
- Removable armrest and head rest covers
- Ergonomic handset with large backlit buttons
- USB port (type A) on handset
- Side pockets on both sides of chair
- Quality German Okin electric actuators
- Battery back up
- Adjustable feet
- Ø38mm (1.5") castors

SPECIFICATIONS

(All specifications may vary (+/- 10mm) as each chair is handcrafted)

Overall (Sitting): 910mm(W) x 1050mm(D) x 1170mm(H)

Seat: 530mm(W) x 520mm(D) x 490mm(H)

Armrest (from top of armrest to seat): 160mm(H)

Backrest: 770mm(H)

Legrest: 590mm (W) x 450mm (L) x 90-210° Angle range

Required distance from wall: 710mm

Maximum User Weight: 158kg

Product Weight:

Base - 54kg / Backrest - 15kg / Combined - 69kg

Please contact **REDGUM Brand** if you require further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

ASK FOR THE BRAND

redgumbrandmobility.com.au