



**10 YEAR WARRANTY**  
terms and conditions apply

**RED UL8-W**

# VENICE



**AIRLINE TRAVEL READY**



## FEATURES

- Powerchair pulls apart for storage and transport
- Mid back contoured seat
- Finger control operation
- High / Low speed settings
- Mid wheel drive
- Battery power indicator / Horn
- Height / angle adjustable seat
- Lift up & adjustable width / height arm rests
- Fold up adjustable height / depth footrest
- Ideal for indoor and outdoor use on even ground
- Airline travel ready - MSDS online
- Lithium powered

## SPECIFICATIONS

- Overall: **725mm (L) x 545mm (W) x 1180mm (H)**
- Turning Circle: **755mm**
- Speed: **6km/h**
- Maximum range: **Up to 12km**
- Safe climb angle: **10°**
- Maximum climb angle: **10°**
- Motor: **2 x 0.5Hp high torque motors (2 Pole)**
- Batteries: **10 Amp minimum (Lithium Iron Phosphate)**
- Charger: **2 Amp Lithium**
- Tyres (Flat free): **Ø149 x 50mm (Front)**  
**Ø2.80 / 2.50 - 4 (Rear)**
- Seat with Controller Weight: **16kg**
- Rear Chassis with Drive Train Weight: **14.5kg**
- Front Chassis with Battery Weight: **13kg**
- Product Weight: **43.5kg - with batteries**
- Maximum User Weight: **113kg - on flat ground**



0° - 113kg  
6° - 100kg  
10° - 90kg

Please contact **SHOPRIDER** for more information on this product

**P: 08 9248 4180**  
**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**

**ASK FOR THE BRAND**

**[shoprider.com.au](http://shoprider.com.au)**

**IMPORTANT NOTE:**  
It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

**SHOPRIDER** strives to provide the highest quality products possible that are subject to change without notice.