



<28.8 km
Range per Charge



170 kg
Weight Capacity



109 cm
Turning Radius



**ZERO
TURN 10**

Effortless Manoeuvrability Equals a Seamless, Comfortable Ride

With contemporary styling and cutting-edge technology, the Zero Turn 10 is the perfect combination of aggressive performance and manoeuvrability. It's EZ Turn Technology and dual motors, allow the scooter to turn on its own axis, making tight turns, narrow halls and doorways a breeze. Combined with Comfort-Trac Suspension (CTS), this scooter maintains a comfortable, stable ride indoors and outdoors.



ZERO TURNIO



Tight Turning Radius



High Back Pillow Top Seat



All Around CTS Suspension

RANGE PER CHARGE (MAX) ^{1,3}	Up to 28.8 km
MAXIMUM WEIGHT CAPACITY	170 kg
TURNING RADIUS ²	109 cm
MAXIMUM SPEED ¹	10 km/h
LENGTH & WIDTH ^{2,5}	L: 122 cm W: 61.6 cm
SEAT TYPE WIDTH & DEPTH (STANDARD)	Pillow top deluxe, H-B w/ sliders W: 49.5 cm (19.5") D: 44.5 cm (17.5")
SEAT-TO-DECK HEIGHT ²	47.5 cm - 52.5 cm
SEAT-TO-GROUND HEIGHT ²	58.4 cm - 63.5 cm
GROUND CLEARANCE ^{2,4}	5 cm at motor

OVERALL WEIGHT (NOT INC BATTERY) ^{2,6}	82 kg
WEIGHT OF HEAVIEST PIECE ²	34 kg (rear section)
SUSPENSION	Front and rear
3W/4W OPTION	4W
FRONT TYRES	9" x 2.75" solid
REAR TYRES	10.75" x 6" solid
BATTERY REQUIREMENTS ⁷	(2) 40AH
BATTERY CHARGER	Off-board, 3.5A
BATTERY WEIGHT ⁸	14 kg
COLOUR OPTIONS	■ Candy Apple Red ■ Black

1. Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
3. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from the theoretical calculations based on battery specification and drive system performance. Test conducted at maximum weight capacity.
4. Does not account for anti-tip clearance.
5. Width of base (not including seat).
6. Includes standard seat weight.
7. Deep cycle (AGM or gel cell type recommended) excluding Lithium-ion models.
8. Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.